

An illustration of a blood smear. It features a large, light pink, irregular shape representing a blood drop. Inside this shape, there are three red blood cells depicted as red, biconcave discs. One cell is at the top left, one is in the middle left, and one is at the bottom center. The background is white.

AI Microscope

-Classify Blood Cells-

Frouke Albrecht



01
Introduction

02
Data

03
Approach

04
Results

05
Conclusions

06
Outlook

The background of the slide features a light pink, organic shape with concentric, wavy lines. Three red blood cells, depicted as red biconcave discs, are scattered within the pink area: one near the top left, one in the middle left, and one at the bottom center.

01

Introduction


Background

- White blood cells important for the **immune system** and defend the body **against infectious disease** and foreign materials
- **5 main types** of white blood cells
- Share commonalities but are **distinct in form and function**



Objective

Automize Classification of
different types of **white**
blood cells in microscope
images



Motivation

- Support medical experts with visual diagnosis
- Increase diagnostic accuracy




The background of the slide features a series of concentric, hand-drawn style circles in various shades of pink. Three red blood cells, depicted as biconcave discs with a dark red color and a lighter red center, are scattered across the left side of the image. The text '02' is positioned on the right side of the pink area.

02


Data

Data




5 main types of white blood cells

Data




5 main types of white blood cells

Data




5 main types of white blood cells

Data




5 main types of white blood cells

Data




5 main types of white blood cells


Examples




Eosinophil



Neutrophil




Lymphocyte



Monocyte

Data

Images



348 images total



Images have been
augmented to about
2500 for each
category




03

Approach





Convolutional Neural Network







input


Feature Importance



1. ConvLayer



2. ConvLayer



3. ConvLayer




04

Results

Results

Accuracy 85%






05

Conclusions

Conclusions



- Developed automated tool for blood cell classification from images
- Classified 4 types of cells
- 85% accuracy
- First step to a bigger problem:
In practice important to detect and count white blood cells



06

Outlook

Outlook: Object Detection



An increased or decreased number of leukocytes indicates the presence of a disorder.

Thank you.



falbrechtg@gmail.com



<https://github.com/froukje/>


CREDITS: This presentation template was created by Slidesgo, including icons by Flaticon, and infographics & images by Freepik.

Data


| Category | |
|------------|---|
| Eosinophil | Responsible for combating multicellular parasites and certain infections; control mechanisms associated with allergy and asthma |
| Neutrophil | Kill bacterias. Migrate toward sites of infection or inflammation |
| Lymphocyte | Lysis of virally infected cells and tumour cells; Release cytokines and growth factors that regulate other immune cells; Immunoregulation and cytotoxicity; Secretion of antibodies |
| Monocyte | Replenishing resident macrophages under normal conditions; migration in response to inflammation signals; differentiation into macrophages or dendritic cells to effect an immune response. |

Data

| Category | Type | Size | % |
|------------|---|---------------|----------|
| Eosinophil | granulocyte/ myeloid | 12-17 μ m | 1 - 3% |
| Neutrophil | granulocyte/ myeloid | 12-15 μ m | 60 - 70% |
| Lymphocyte | mononuclear leukocyte (agranulocyte)/ lymphoid | 6-10 μ m | 25 - 40% |
| Monocyte | mononuclear leukocyte (agranulocyte)/ myeloid | 12-10 μ m | 2 - 10% |




Principal Component Analysis




Principal Component Analysis


Explained Variance



Eigenvectors



Principal Component Analysis



Convolutional Neural Network

